

129

750

TK 51.285

KFKI-75-54

I. MANNO

COMPUTER CONTROL OF BIS SPECTROMETER

*Hungarian Academy of Sciences*

CENTRAL  
RESEARCH  
INSTITUTE FOR  
PHYSICS

BUDAPEST



1975 NOV 21

1975 11/21

2017

THE UNIVERSITY OF CHICAGO PRESS  
50 EAST LEXINGTON AVENUE  
NEW YORK, N.Y. 10017  
TEL: 212 850 6645  
WWW.CHICAGO.PRESS.COM

KFKI-75-54

COMPUTER CONTROL OF BIS SPECTROMETER

I. Manno

Central Research Institute for Physics, Budapest, Hungary  
High Energy Physics Department

ISBN 963 371 060 X

## ABSTRACT

This paper has the purpose first to familiarize the interested reader with the computer control of the BIS spectrometer, and then to help him to use the Control Program.

## Аннотация

Контрольная система была выработана для спектрометра БИС-а, работающего на линии с ЭВМ ТРА-1.

## KIVONAT

Ez a leírás a BIS spektrométer számítógépes ellenőrzésével ismerteti meg az olvasót, és egyben segítséget nyújt a felhasználóknak a Control Program használatához.

## CONTENT

I.	Introduction . . . . .	1
II.	BIS Spectrometer . . . . .	1
III.	Control Equipment. . . . .	2
IV.	Control Program. . . . .	2
V.	Commands . . . . .	4
VI.	FOCAL. . . . .	5
VII.	Input/Output . . . . .	5
VIII.	Diagnostic Messages. . . . .	5
IX.	Literature . . . . .	6
X.	Appendices	
	A. ODT Command Summary. . . . .	7
	B. FOCAL Command Summary. . . . .	8
	C. Loading Procedure. . . . .	10



## 1. INTRODUCTION


The computer control of the BIS spectrometer allows to recognize and to repair the defects of the spectrometer easily. Such a way expensive proton synchrotron time is saved and statistics of events may be collected for experiments in a minimal time.

The described control system allows to change the BIS spectrometer from one experiment to the other easily.

## II. BIS SPECTROMETER

Several important experiments were carried out with the BIS spectrometer at the 76 GeV proton synchrotron of the high energy laboratory ИФВЭ\* in Serpukhov. The most important experiments are the measurement of  $K^0$ -regeneration, the investigation of  $K_S \rightarrow \mu\mu$  decay and the search for charmed particles.

The BIS spectrometer (Fig.1.) consists of magnets, spark chambers, proportional chambers, scintillator counters,  $\mu$ -detectors and electron-detectors.



\* Институт Физики Высоких Энергий

The spectrometer works on line to a BESM-3 computer. This computer records the events onto magnetic tapes. These recorded tapes are evaluated off-line on computers: BESM-6, CDC-6200 etc.

### III. CONTROL EQUIPMENT

The BIS spectrometer is on line to a TPAi small computer too /Fig.2./




Fig.2.

The installed TPAi computer consists of a 16K core memory, 32K NC-245i Disc Unit, PERFORMOM 30 punch, FS 1501 reader, ASR33 Teletype and NE-601/i TV-display.

There is an interface with 256 12-bit word buffer memory between the spectrometer and the TPAi small computer, Less than 0.3 msec is needed to fill up the buffer memory. The interface drives a point display too.

During the measurement the data are transferred from the BIS spectrometer via the buffer memory into the computer core memory. The

Control Program selects, controls and evaluates the measured data according to the commands given. There are commands to type or to display the results.

### IV. CONTROL PROGRAM

The Control Program consists of program segments. The program segments are called by commands typed at the Teletype in the form of maximum four characters. All commands are executed by typing the RETURN key or a comma directly after the commands /Fig.3./.

The Control Program is split into two sections: initialization and execution.

a.) During initialization the operator may change the desired parameters of the Control Program.

Using the incorporated ODT debugging program the programmer can run his binary subroutine on the computer, control its execution, examine






Fig. 3.




Fig. 5.




Fig. 4.

registers, change their contents and make alterations to his program by typing at the Teletype keyboard. The programmer may have the corrected program in binary. Finally he can incorporate his program into the Control Program. /The precise capabilities and commands of ODT are detailed in [2], [6] and Appendix A./

The Control Program may type its command set /Fig. 6 /.

b./ During execution the operator may direct the evaluation of the data by typing commands detailed below.

NOTE: The measured data overlay the command table and the ODT.

v. COMMANDS

*CMND		KILL	Clear the data buffer. The buffer contains the histograms and the views of the events.
COMMANDS:			
KILL	CLEAR BUFFER		
NCH	CHAMBER NO?	NCH	The Control Program waits for the serial number of the desired spark chamber. The operator must type the serial number in floating point format.
NEV	EVENT NO?		
MRK	MARKER?		
X	X VIEW		
Y	Y VIEW		
LHT	DISPLAY HISTOGRAM		
LCH	DISPLAY CHAMBERS		
S	STOP		
C	CONTINUE		
CLR	CLEAR DISPLAY		
SLT	SELECT		
FCL	START FOCAL		
GO	START RUN		
CMND	COMMANDS	NEV	The Control Program waits for the serial number of the desired event in a accelerator cycle. The operator must type the serial number in floating point format.
ODT2	START ODT		

\*

Fig. 6.

MRK	The Control Program waits for the marker. The operator must type the place of the marker in floating point format.
X	X view.
Y	Y view.
DHT	Display histogram according to the NCH, NEV, MRK, X and Y commands.
DCH	Display the view of the event according to the NEV, X and Y commands.
S	Stop execution.
C	Continue execution.
CLR	Clear the display buffer.
STL	Select program to evaluate data. The user may write no more than 7 evaluating programs.
FCL	Transfer control to the FOCAL on line interpreter.
CMND	Type the command set /Fig. 6./.
ODT2	Transfer control to the ODT debugging program.

## VI. FOCAL

The Control Program contains the FOCAL on line interpreter./You may find detailed description of FOCAL in [3], [5], [10], and Appendix B./

Using FNEW function [7] FOCAL programs may read data from the core memory or may store data there in the form of one word fixpoint format.

The call of FNEW is as follows:

```
SET Z = FNEW /WR, F, LC,V/
```

Parameters:

$$WR = \begin{cases} 1 & \text{reading} \\ -1 & \text{storing} \end{cases}$$

F = The serial number of the field /4K/

LC = Decimal address of the word in field F

V = Value to store in case of writing

Any parameter may be a number, a variable or an expression.

The user may write programs in FOCAL to evaluate or to display [10] measured data.

The FCL command transfers the control to the FOCAL.

## VII. INPUT/OUTPUT

Information is transferred between peripheral devices and the TPA1 small computer through program interrupt. The devices signal the computer when they are ready to transfer information, the program will then interrupt its normal flow and jump to a routine /SERV/ to process information, after which it will return to the point in the main program at which it was interrupted. Thus the transfers are device-initiated but under program control.

## VIII. DIAGNOSTIC MESSAGES

Error messages of FOCAL are described in [3] and [5].

When an error occurs during program execution, error message is typed on the Teletype. The Control Program types the word ERROR AT and the location where the error is encountered /Fig.7./.




Fig. 7.

IX. LITERATURE

- [1] PDP8/e 1 PDP8/m Small Computer Handbook 1972. Digital Equipment Corporation, 1972, Massachusetts, USA
- [2] Introduction to Programming, DEC, 1970
- [3] Programming Languages, DEC, 1970
- [4] Edusystem Handbook, DEC, 1973
- [5] FOKAL, KFKI TPA-IY-01-MA, TPA-IY-02-MA, TPA-IY-03-MA
- [6] ODT, KFKI 5114
- [7] I.Manno: FOCAL in On-line Data Processing, 1975, DECUCOPE, Massachusetts
- [8] I.Manno: SOFT, a Supervisor of both FOCAL and On-line Measuring Terminals, 1974, KFKI-74-53
- [9] OS/8 Software Support Manual, DEC, 1973
- [10] J.Eszenszki and I.Manno: Display-FOCAL, 1975, KFKI

APPENDIX A

ODT COMMAND SUMMARY

nnnn/           Open location designated by the octal number nnnn.  
/                Reopen last opened location.  
RETURN key      Close previously opened location.  
LINE FEED key   Close previously opened location and open the next sequential  
                  one for modification.  
↑ (SHIFT/N)     Close location, take contents of that location as a memory  
                  reference and open it.  
+ (SHIFT/O)     Close location, open indirectly.  
Illegal character Current line typed by user is ignored, ODT, types  
                  ?(CR/LF).  
nnnnG           Transfer program control to location nnnn.  
nnnnB           Establish a breakpoint at location nnnn.  
B                Remove the breakpoint.  
A                Open the location in which the contents of AC were stored when  
                  the breakpoint was encountered for modification.  
LINE FEED key   Open the location in which the contents of L were stored when  
                  the breakpoint was encountered for modification.  
C                Proceed from a breakpoint.  
nnnnC           Continue from breakpoint and iterate past the breakpoint nnnn  
                  times before interrupting the user's program at the breakpoint  
                  location.  
M                Open the search mask.  
LINE FEED key   Open the lower search limit.  
LINE FEED key   Open the upper search limit.  
nnnnW           Search one portion of core as defined by the upper and lower  
                  limits for the octal value nnnn using the mask.  
T                Punch leader.  
nnnn;mmmmP      Punch a binary core image defined by the limits nnnn and mmmmm.  
E                Punch checksum and trailer.

APPENDIX B

FOCAL COMMAND SUMMARY

Command	Abbreviation	Example of Form	Explanation
ASK	A	ASK X,Y,Z	FOCAL types a colon for each variable, the user types a value to define each variable.
COMMENT	C	COMMENT	If a line begins with the letter C, the remainder of the line will be ignored.
CONTINUE	C	CONTINUE	Dummy line
DO	D	DO 4.1	Execute line 4.1; return to command following DO.
		DO 4.0	Execute all the lines beginning with 4. Return to the command following DO also when a RETURN is encountered
		DO ALL	Execute the whole program
ERASE	E	ERASE	Erase the symbol table
		ERASE 2.1	Delete line 2.1
		ERASE 2.0	Erase all the beginning with 2.
		ERASE ALL	Delete all user input
FOR	F	FOR I=X,Y,Z; /commands/ FOR I=X,Z; /commands/	Execute the commands following unless I is greater than Z. The initial value of I is X and at the end of each cycle Y is added to I.
GO	G	G	Start indirect program at the lowest line number
GO?	G?	GO?	Start at the lowest line number and trace entire indirect program until another? is encountered or until the completion of program.
GOTO	G	GOTO 3.4	Start indirect program /transfer control to line 3.4/ Must contain an argument.
IF	I	IF/X/Ln,Ln,Ln IF/X/Ln,Ln /commands/ IF/X/Ln /commands/	Where X is a defined identifier a value, or an expression, followed by one to three line numbers. If X is less than zero control is transferred to the first line number, if X is equal to zero to the second line number, if X is greater than zero to the third line number.

MODIFY	M	MODIFY 1.15	Enables editing of any character on line 1.15. The (CTRL/L), (CTRL/BELL), RUBOUT, + , LINE FEED and RETURN keys may be used,
QUIT	Q	QUIT	Return control to the user
RETURN	R	RETURN	Terminates DO subroutines, returns to the original sequence
SET	S	S A=5/B-C	Defines identifiers in the symbol table
TYPE	T	TYPE A+B-C	Evaluate expression and type out = and the result in current output format.
		TYPE A-B,C/D	Compute and type each expression separated by commas.
		TYPE "TEXT"	Types texts. May be followed by ! to generate new line, or # to generate carriage return.
WRITE	W	WRITE	FOCAL types out the entire indirect program.
		WRITE ALL	
		WRITE 1.0	FOCAL types out all the lines beginning with 1.
		WRITE 1.1	FOCAL types out line 1.1

OPTIONS:

OPTION	R	High-speed reader input.
OPTION	K	Keyboard input.
OPTION	P	High-speed punch output.
OPTION	I	Interpretive Input mode and Numeric Output mode.
OPTION	C	Character Input/Output mode.
OPTION	X	Suppress the colon printout.
OPTION	:	Restore the colon printout.
OPTION	S n	Set the number following "s" as the character code for an extra input terminator. Decimal value ASCII code.
OPTION	S	Delete the special terminator.
OPTION	N	No echo.
OPTION	E	Restore the echo.
OPTION	M	Start the Disk Monitor.

APPENDIX C

LOADING PROCEDURE

There are two ways to load the Control Program.

a./ The Control Program is punched on two paper tapes in binary format. To load the paper tapes the FIELD LOADER must be in core.

Place binary tape 1 into the reader, set the switches to 7777 then press the LOAD ADDRESS key, set the switches to 3777 then press the CLEAR and CONTINUE keys in that order. When the tape has been read, the status of the accumulator lights will indicate any error in loading. If the lights are out, the loading was successful; if any lights are on, there was a checksum error and the tape must be reread. Repeat the described procedure with binary tape 2.

The Control Program starts at 200 in field 2.

b./ The user may load the Control Program using the Disc Monitor System. This procedure is shown on Figure 8.



.PIP  
\*OPT-L

\*IN-S:

FE=0033

NAME TYPE ELK

AF  
PALD.SYS (0) 0037  
PIP .SYS (0) 0025  
LOAD.SYS (0) 0011  
.CL .SYS (0) 0007  
ELIT.SYS (0) 0015  
ASC .ASCII 0103  
ODT .USER(2) 0004  
FLTG.USER(2) 0011  
FCL .USER(2) 0037  
OL1 .USER(2) 0004  
OL2 .USER(2) 0016

\*OPT-

.CALL OL1  
.CALL OL2  
.CALL OLT  
.CALL FLTG  
.CALL FCL  
\*GMND

COMMANDS:

KILL CLEAR BUFFER  
NCH CHAMBER NO?  
NEV EVENT NO?  
MEK MARKER?  
X X VIEW  
Y Y VIEW  
DHT DISPLAY HISTOGRAM  
LCH DISPLAY CHAMBERS  
S STOP  
C CONTINUE  
CLR CLEAR DISPLAY  
SLT SELECT  
FCL START FOCAL  
GO START RUN  
GMND COMMANDS  
OIT2 START OLT

\*GO  
\*FCL

C-FOKAL, 1971


APPENDIX D

EXAMPLE


BEAM DISTRIBUTION AT THE 1ST CHAMBER Y VIEW  
MEASURED BY THE BIS SPECTROMETER  
DATE: 16. V. 1975. SERPUCHOV

BIN CONT ENTRIES= .514000E+04 MEAN= .309813E+02

1	1
2	0
3	2
4	1
5	1
6	2
7	5 X
8	5 X
9	5 X
10	14 XXX
11	24 XXXXXX
12	23 XXXXXX
13	29 XXXXXXXX
14	32 XXXXXXXX
15	43 XXXXXXXXXXX
16	38 XXXXXXXXXXX
17	59 XXXXXXXXXXXXXXXX
18	79 XXXXXXXXXXXXXXXXXXXX
19	91 XXXXXXXXXXXXXXXXXXXXXXXX
20	106 XXXXXXXXXXXXXXXXXXXXXXXXXXXX
21	130 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
22	139 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
23	160 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
24	171 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
25	187 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
26	193 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
27	216 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
28	227 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
29	231 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
30	239 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
31	240 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
32	237 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
33	238 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
34	219 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
35	218 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
36	199 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
37	185 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
38	176 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
39	143 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
40	135 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
41	128 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
42	107 XXXXXXXXXXXXXXXXXXXXXXXX
43	77 XXXXXXXXXXXXXXXXXXXXXXXX
44	86 XXXXXXXXXXXXXXXXXXXXXXXX
45	73 XXXXXXXXXXXXXXXXXXXXXXXX
46	60 XXXXXXXXXXXXXXXXXXXXXXXX
47	44 XXXXXXXXXXXXXXXX
48	28 XXXXXXXX
49	31 XXXXXXXX
50	22 XXXXX
51	9 XX
52	4 X
53	6 X
54	9 XX
55	5 X
56	4 X
57	3



62244



Kiadja: a Központi Fizikai Kutató Intézet  
Felelős kiadó: Pintér György, a KFKI  
Részecske- és Magfizikai Tudományos  
Tanácsának szekcióelnöke  
Szakmai lektor: Telbisz Ferenc  
Nyelvi lektor: Sebestyén Ákos  
Példányszám: 450 Törzsszám: 75-972  
Készült a KFKI sokszorosító üzemében  
Budapest, 1975. augusztus hó

